

The **Overnighter™** improves upon the comfortable performance of the Hickory Springs ClassicRest™ mechanism. The Overnighter™ provides a number of additional features which make it the ideal mechanism for residential use. And like all Hickory Springs sleeper mechanisms, Overnighter™ units are performance-tested thousands of times to ensure a superior level of dependability and customer satisfaction.

- One-piece crossover leg tube - provides extra durability and stability when the unit is open and is contoured to protect the consumer when opening.
- Semi-tubular heavy gauge steel perimeter - provides added durability.
- Galvanized helicals - provide additional support and durability to polypropylene seating and sleeping deck.
- Headlock™ locking TV headrest system - prevents bed from being closed with raised headrest.
- Synchronized pivotal action hinge - (newly redesigned) allows for smoother operation and eliminates sliding parts which could damage bedding. Also allows for maximum clearance of upholstered arms while folding or unfolding.
- Nylon washers in critical rivet joints - reduce wear and ensure easy operation.
- Smoothed exposed edges - prevent upholstery damage or consumer injury.
- One-piece front leg system - ensures a consistent seat pitch from side to side.
- Positive locking leg tube - allows the unit to be opened easily from either side of the mechanism.
- Contoured "saddle seat" in seat section - provides added consumer comfort.
- Sleeper cavity - accommodates either innerspring or foam mattress.
- Tilt-forward feature - allows for easy cleaning under the sofa.
- Charcoal (black) color - coordinates with most styles and colors of upholstered fabrics and provides a shadowing effect between cushions.

Standard OverNighter Model 900

- BedSecure™ feature - prevents the bed from tilting or toppling at the head section, reducing the chance of consumer injury
- Houses up to 5" innerspring mattress 75½" long.